

Shanghai Tenry Pharmaceutical Co., Ltd


No.151, Lane.1236, Xinsiping Road, Situan Town, Fengxian District,
Shanghai, 201413, P.R.China

It is certified that the manufacturer's technical file and the PPE product detailed on
page 2 have been assessed and found to be in accordance with

Regulation (EU) 2016/425 Module B, EU type-examination

This certificate is valid from 27 October 2020 until 27 October 2025
1. Certified since 27 October 2020

Authorised by




FINAS
Finnish Accreditation Service
S003 (EN ISO/IEC 17065)

SGS FIMKO OY, Notified Body 0598

Takomotie 8, FI-00380, Helsinki, Finland
t +358 9 696 361 f +358 9 692 5474 www.sgs.com

Page 1 of 2



Shanghai Tenry Pharmaceutical Co., Ltd

Regulation (EU) 2016/425 Module B, EU type-examination

Issue 1

PPE Product

TENRY (logo) TR-ZD01 particle filtering half mask, consisting of a white five layer (polypropylene/ polypropylene/ polypropylene / polypropylene / polypropylene) disposable face mask, with nose clip, nylon/spandex ear loop and adjustable hook.

It is certified that the manufacturer's technical file and the above mentioned PPE have been assessed and found to meet the applicable Essential Health and Safety Requirements in Annex II of Regulation (EU) 2016/425 Personal Protective Equipment

The following have been applied:

EN 149:2001+A1:2009 (Respiratory protective devices - filtering half masks to protect against particles) device classification: FFP2 NR.

This certificate is issued on the strict condition that appropriate checks on manufactured PPE, as detailed in Article 19 (c) of the Regulation are implemented and maintained while the model is in production

Certification is based on technical file reference:
Particulate Respirator Protective Mask /TR-ZD01,
version 1, dated: 2020-09-17.

SGS Reference Number UK/CRS 242225.

This certificate remains the property of SGS Fimko Oy to whom it must be returned on request